



RESTORING PET MOBILITY:

21 BENEFITS OF WIMBA PET ORTHOTICS




About WIMBA

Step into a world where your pet's comfort and mobility come first. Explore the main 21 benefits of WIMBA orthotics, meticulously crafted to cater to your pet's unique needs. With advanced materials and cutting-edge technology, WIMBA orthotics ensure lasting support and comfort, redefining how pets thrive through every step of their journey.

We specialize in creating custom orthotics and prosthetics solutions for pets, aiming to enhance their mobility and overall comfort. Our orthotics are meticulously crafted using advanced materials like nylon and TPU, renowned for their durability and flexibility. Designed to be waterproof and lightweight, WIMBA orthotics ensure pets can maintain their active lifestyles without discomfort.

Whether you're a pet parent looking to enhance your companion's quality of life or a veterinary professional interested in offering advanced orthotic solutions, WIMBA is here to support you. Join us in redefining pet care with our innovative orthotic solutions.



21 Benefits of WIMBA Pet Orthotics

Lightweight



4x lighter construction, made with ultralight printing materials like nylon and TPU, ensures that pets can move freely and comfortably, minimizing additional strain or discomfort.



Durable



WIMBA orthotics are crafted from high-quality, durable materials that withstand daily wear and tear, ensuring long-lasting support for pets when used properly.



Custom-fit



Utilizing HP MultiJet Fusion 3D printing and AI technologies, WIMBA orthotics are custom-fitted with pinpoint accuracy, providing optimal support and enhancing the effectiveness of the treatment.



21 Benefits of WIMBA Pet Orthotics



Precise Scanning



Efficient scanning with the help of WimbaSCAN, accurate to 0.2 mm, and manufacturing processes reduce fitting time and improve accuracy, enhancing overall patient care.

Waterproof



Designed to be waterproof and easy-to-clean for optimal hygiene, WIMBA orthotics let pets enjoy their activities without compromising effectiveness and can be used during water treadmill rehabilitation.



Adaptive Rehabilitation




The unique WimbaTOKENS system allows for customizable resistance and range of motion adjustments, ensuring that the orthotic device adapts to the pet's rehabilitation progress and specific needs.

21 Benefits of WIMBA Pet Orthotics

Cushioned Inserts



Soft and breathable, 3D-printed elastic inserts provide additional comfort and support, helping to prevent sores, ensure a snug fit for the pet, and ensure airflow around surgical wounds.



Affordable



WIMBA orthotics, especially the stifle brace, can be an affordable option for pet parents with financial constraints that currently limit access to surgery.



Surgical Alternative



WIMBA orthotics, especially the stifle brace, can be a long-term support and alternative for pets who aren't candidates for surgery.



21 Benefits of WIMBA Pet Orthotics



Preventive



After TPLO or TTA surgery, dogs face a 22%-54% risk of ACL rupture in the opposite leg within 6 to 17 months, and the WIMBA Stifle orthotic serves as a preventive solution.

Mess-Free



WIMBA replaces messy plaster casting with digital scanning and 3D printing, streamlining the fitting process for pets and vets, taking only 15-20 min to scan, customize, and order via the WimbaAPP.



Enhanced Mobility



Designed to restore and improve mobility, WIMBA orthotics help pets regain their movement capabilities and enhance their quality of life.



21 Benefits of WIMBA Pet Orthotics

Eco-Friendly



WIMBA's use of sustainable materials and precision 3D printing technology reduces waste, supporting eco-friendly practices while delivering top-quality orthotic solutions.



Modern Design



With a sleek and modern design, WIMBA orthotics not only function well but also look appealing, making them a preferred choice among pet parents and vet professionals.



Secure Fit



Adjustable straps and secure closures ensure that WIMBA orthotics stay in place, providing consistent support throughout the day.



21 Benefits of WIMBA Pet Orthotics



Rapid Turnaround



With most WIMBA orthotics delivered in 2 weeks or less, pets receive timely treatment, reducing wait times and accelerating recovery.

Comprehensive Care



WIMBA orthotics address a wide range of orthopedic issues, providing more holistic and effective care to patients.



Proven Effectiveness



Backed by veterinary professionals, WIMBA orthotics have a track record of improving mobility and recovery in pets.



21 Benefits of WIMBA Pet Orthotics

Optimal Comfort



The ergonomic design and soft padding of WIMBA orthotics ensure that pets remain comfortable during use, reducing the risk of discomfort or irritation.



Ease of Use and Adjustment



WIMBA orthotics are designed for easy application and removal, while the design allows for easy adjustments with interchangeable inserts, tokens, and straps to adapt to the pet's changing needs.



Customer Support



Access to expert advice and WIMBA support ensures that veterinarians can effectively use and fit WIMBA orthotics for their patients.



What's Next?

Are you a Pet Parent?

Ready to give your pet the gift of enhanced mobility and comfort? Contact us today to order WIMBA orthotics or check our map of current WIMBA providers to locate one near you and order directly.

[CONTACT US](#)



Are you a Vet Professional?

Interested in offering plaster-free, custom-made orthotic solutions in your practice? Learn how to become a WIMBA provider and start providing advanced care with WIMBA orthotics. Contact us now to get started on your journey to becoming a WIMBA provider and enhancing your practice's orthotic offerings.

[BECOME A PROVIDER](#)

